# Abzugsschrank ASEX420

## Extraction cabinet ASEX420

### Sorbonne ASEX420

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### Instruction for use

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The safety of persons and equipment in potentially explosive environments depends on compliance with all safety regulations. The prerequisite for the assembly, maintenance, repair work on the spark resistant equipment is therefor the exact knowledge of the applicable regulation and provision in particular relating to

- The general recognized rule of technology
- The national accident prevention and safety regulations
- The regulation and provision applicable for the spark resistant equipment
- The information in this manual
- The declaration on the identification plate and the information signs

### **Dear Customer**

Thank you for choosing this high-quality product from BDT.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

These written instructions cannot cover all possible situations that can show up during use. For this reason, our hotline is available to answer all questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, to study the enclosed current instructions for use or to refer to our website at [www.bdt-dental.de](http://www.bdt-dental.de).

### **Manufacturer**

BDT GmbH • Industriestraße 27 • 77656 Offenburg • Germany  
Telephone +49 781/58688  
Fax +49 781/58664  
E-mail [info@bdt-dental.de](mailto:info@bdt-dental.de)

### **With reference to our quality standards**

BDT ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

---

# Instruction for use Extraction cabinet ASEX420

Article no. 10046

## 1. Safety information

### 1. Correct usage

- The extraction cabinet is designed solely for the storage and processing of hazardous materials/ substances in the dental / orthodontic laboratory. Any other or further use beyond these parameters is considered improper use.
- The correct use of the unit includes compliance with all instructions in the instruction manual.
- The extraction cabinet may only be operated in an undamaged condition
- Structural changes to the extraction cabinet can lead to hazards, impairments, loss of spark resistant protection and in any case lead to loss of declared conformity
- The spark resistant components may only be replaced by original spare parts
- When a BDT product is finally no longer in use the disposal regulations of that particular country apply. BDT or the dental trade supplier can be contacted to answer questions regarding the correct disposal of any specific product.

### 2. Danger symbols and markings

In the instructions for use, the following symbols are used to indicate sources of danger to the user:

#### Warning

This symbol indicates direct danger to the life or health of the user. Disregarding this warning may involve a serious risk to health.

#### Caution

This symbol indicates important safety instructions. Disregarding this warning may cause damage to the unit and/or other property.

### 3. Disposal

This symbol indicates that products marked with it should not be disposed of together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through BDT or the dental trade.

---

## 2. Technical data

productions series	ASEX420
Certificate no.	IECEX TPS 20.0023X
spark resistant	Ex h IIA T60 °C Gc/-
	 II 3G/- Ex h IIA T60 °C Gc/-
power supply	230 V 50/60 Hz
rated power	200 W
measurements (w x h x d)	1000 x 1290 x 500 mm
weight	65 kg
protection category	I
degree of protection	IP 55
flow volume	420 m <sup>3</sup> /h
pressure	stat (max) 340 Pa
 operational temperature range	10 °C to 40 °C (50 °F to 104 °F)
utilization	only indoors
height above sea level	max. 1500 m above sea level
 air humidity	> 55 %
maximum load power outlet	1380 W

## 3. Standard conformity

The extraction cabinet complies with the listed standards, in the separately enclosed declaration of conformity. References to standards and guidelines in this manual always refer to the current version. Additional supplements (for example, dates) must be observed.

---

## 4. Commissioning of the Extraction cabinet

If, contrary to expectations, you notice a transportation damage to your device, immediately notify BDT or the specialist dealer from whom you obtained the device.

Before commissioning, make sure that the mains voltage complies with the operating voltage specified on the rating plate (rear of the device).

After setting up the extraction cabinet, an exhaust duct leading into the open air must be connected to the extraction cabinet. The diameter of the outlet on the fan is 100mm. The maximum length of the exhaust duct may not exceed 12m and for each 90° bend 1m must be deducted. A suitable exhaust duct with explosion approval must be used.

The position of the vent opening to the outside must be free of ignition sources (e.g., electrical appliances, switches, or hot, glowing surfaces).

Now connect the device to the electrical mains (230V).

The left switch turns on the light installed in the extraction cabinet (Fig. 1).

The right switch turns on the fan installed in the extraction cabinet.

### Operation of the extraction cabinet

The extraction cabinet is designed for a 24-hour operation.

If an electrical device is operated in the extraction cabinet, the supply line of the respective device can be led out through the cable grommet on the left side. You will find an additional output for one device with cold-device plug on the distribution box of the extractor cabinet. There the device can be connected with the supplied cold-device plug extension. The maximum allowed power consumption is 1380 watts.

The electrical equipment inside the fume cupboard must either be explosion approved or constructed so that no sparks can be generated during installation and operation.

If the unit / container located in the extraction cabinet emits a fuel that can generate an explosive atmosphere with air or if such a process cannot be ruled out, the extraction cabinet is in continuous operation for the duration of the aforementioned fuel emission with a follow-up time of 5 minutes operate.

The maximum permissible fuel quantity is 8 litres.

Heat sources that are warmer than 60 ° C or which heat the atmosphere of the extractor cabinet to more than 60 ° C must not be operated.

When using the lighting of the extractor cabinet the fan must be switched on.

The use of lighting without the fan switched on is strictly prohibited.


Fig. 1

## 5. Maintenance / servicing

- The regulations governing the maintenance / servicing of electrical equipment in potentially explosive atmospheres, such as EN / IEC 60079-17, must be adhered to.
- Before opening the housing, make sure that there is no voltage present or take suitable protective measures.
- The cleaning of the extraction cabinet must be cleaned only with a damp cloth due to possible static charge inside and outside.

## 6. Repair / Maintenance / Change

Maintenance work / repairs may only be made with original parts.

In the event of damage to the spark resistant components (motor, switch and light), only replacement with original products. The decision is to be made only by trained, qualified personnel familiar with the relevant hazards or consulting with BDT.

Repairs relating to spark resistant's of the product may only be carried out by a qualified electrician in accordance with national regulations (EN / IEC 60079-19).

Rebuilding or modifications to the fume cupboard are not permitted

## 7. Problemsolving

Problem	Troubleshooting	Repair
Horns off the distribution box X1	Overvoltage protection has tripped	Replace the overvoltage protection F5
Extraction cabinet doesn't switch on the light or fan	Plug not plugged in	Insert the plug in to the socket
	Fuse F3 has tripped	Reset fuse in F3 in distribution box x1
Extraction cabinet switches on the light but not the ventilation	Motor protection switch F4 has tripped	Reset motor protection switch F3 in distribution box X1
Extraction cabinet works, however the connected electrical device to the socket doesn't work	Fuse F2 has tripped	Reset fuse F2 switch in the distribution box
	Connected device is defect	Contact manufacturer of the product
Extraction cabinet doesn't switch on the light or fan and however the connected electrical device to the socket doesn't work	Mainfuse F1 tripped	Reset fuse F1 switch in the distribution box

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