POLYMERISATIONSWANNE

für 4 Küvetten

POLYMERIZATION BATHS

for 4 cuvettes

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Instruction for use

Prüfbericht / test document

DIN EN 61010

Polymerisationswanne / Polymerization bath IVO - 4 Rev2.0

Seriennummer / Serial number.....	_____
Spannung / Voltage setting.(delete).....	230 Volt / Voltage
Verbrauchstrom / Consumption power (delete).....	_____ Ampere / power
Schutzleiterwiderstand 25 A Prüfstrom / Protective conductor resistance 25 A testing current	> _____ Ohm
Isolationswiderstand / Insulation's resistance	< _____ M - Ohm
Dichtigkeit Luft / Consistence air.....	<input checked="" type="checkbox"/> OK
Dichtigkeit Wasserzulauf / Consistence water intake	<input checked="" type="checkbox"/> OK
Dichtigkeit Wasserablauf / Consistence water discharge.....	<input checked="" type="checkbox"/> OK
Dichtigkeit Becken / Consistence water discharge.....	<input checked="" type="checkbox"/> OK
Temperaturregelung / Temperature regulation.....	<input checked="" type="checkbox"/> OK
Trockenlaufschutz / Dry-running protection	<input checked="" type="checkbox"/> OK
Wasserstandregelung / Water level regulation.....	<input checked="" type="checkbox"/> OK

Zubehör / Supplies:

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Dear Customer

Thank you for choosing this high-quality product from BDT.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

These written instructions cannot cover all possible situations that can show up during use. For this reason, our hotline is available to answer all questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, to study the enclosed current instructions for use or to refer to our website at www.bdt-dental.de.

Manufacturer

BDT GmbH • Sonnhalde 1 • 77749 Hohberg • Germany

Telephone +49 7808 914996-0

Fax +49 7808 914996-99

E-Mail info@bdt-dental.de

With reference to our quality standards

BDT ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the operation of the products.

Responsibility for failures cannot be taken, as we have no influence on the operation on site.

BDT cannot be held responsible for products produced using the equipment due poor operation by the user.

Instruction for use

Polymerization baths

Artikelnr. 10003

1. Safety information

1. Correct usage

The polymerization bath is intended only for polymerizing plastic prostheses in the dental laboratory.

Any other use shall be regarded as not intended. Intended use shall also include compliance with any and all information contained in the instructions for use.

The corresponding regulations of the respective country must be complied with once the BDT equipment is finally taken out of service. BDT or the specialized dental retailer can answer any questions with regard to proper disposal of the equipment.

2. Danger symbols and markings

In the instructions for use, the following symbols are used to indicate sources of danger to the user:

Warning This symbol indicates direct danger to the life or health of the user. Disregarding this warning may involve a serious risk to health.

Warning This symbol indicates that the device produces hot surfaces and burn hazard.

Caution This symbol indicates important safety instructions. Disregarding this warning may cause damage to the unit and/or other property.

3. Disposal

This symbol indicates that products marked with it should not be disposed of together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through BDT or the dental trade.

2. Description and usage

The polymerization bath is used for polymerizing plastic prostheses in the dental laboratory.

The polymerization bath consists of a stainless steel bath with electronic water level control. The temperature can be continuously controlled. As a result of indirect heating there is no furring-up of scaling.

The SR pressure attachments can be connected directly to the polymerization bath by means of the four air plug sockets.

Time measurement is carried out by means of the four adjustable timers.

NOTE: Please note that optimal polymerization can only be ensured with regularly serviced equipment. In this case national regulations must also be taken into consideration.

3. Technical data / environmental conditions

Mains connection	220	V	50/60Hz
Rated power	2500	W	
Protection class	1		
Degree of protection	IP33		
Fuse	6.3 x 32 mm / 12.5 A / 250 V / slow-blow		
Air pressure	NW7	pressure max.	10 Bar
Hose connection	3/4"		
Discharge hose	Ø22mm		
Dimensions (WxHxD)	530x300x470mm		
Weight	14kg		

Operable temperature range	10°C up to 40°C (50°F up to 104°F)
Use	indoor
Height above sea level	max. 1500 above sea level
Air humidity	10% - 80%, non-condensing

4. Commissioning

Immediately notify BDT or the specialized dealer where you purchased the equipment, should the equipment appear to be damaged in transit.

Prior to commissioning make certain that the mains voltage agrees with the operating voltage indicated on the model identification plate (back of the equipment).

After setting up the equipment (distance to the wall and other equipment at least 10 cm) the water inlet hose must be attached to the water supply 3/4" (max 10 bar)

Afterwards the waste water hose must be attached to the waste system. Now the air plug is inserted into the compressed air supply (max. 10 Bar). After the equipment is connected to the electrical mains (230 V plug with earth conductor) the standpipe is inserted into the overflow valve. The equipment is switched on by pressing the main switch. The water supply is now opened and the basin automatically fills up to the adjusted water level. Heating is only possible when the prescribed water level has been reached. The required temperature which must be **below the boiling point** is now adjusted on the thermostat. To avoid the development of steam, the Floating ball can be used. (picture 1)

For safety reasons the heating switches off automatically if the water level is too low.

Distance to the wall and other equipment at least 10 cm, so that a heat accumulation can be avoided.


The plug is the disconnecting device and must be accessible

The water intake must be turned off each evening and the equipment disconnected from the electrical mains.

Do not fill in combustible liquids (fire hazard).

picture 1

5. Care and maintenance

If necessary wipe off the housing with a dry or dampened cloth. Do not use any strong solvents.

The stainless steel bath should preferably be cleaned with a basic stainless steel cleaner and Breglit stainless steel care spray.

The spare parts list is located at the end of the manual or under www.bdt-dental.de.

The water level sensor must be cleaned once per month with a dampened cloth if lime and gypsum residues are cleaned!

Water intake and drip cap blanket has to be tested of pressure tightness.

Repairs should only be carried out by authorized technical personnel.

Only open the housing if the equipment has been disconnected from the electricity mains.

7. Error search

Problem	Possible Cause	Remedy
Green control lamp not lit	Mains cable unplugged	Plug in mains cable
	Fuse defect	Replace fuse
	No mains voltage	Check mains voltage
Water always flows in	Solenoid valve does not close	Clean solenoid valve
	Standpipe is improperly or not installed	Check proper fit of standpipe
	Sensor set at wrong height	Set correct sensor height
	Dirty sensor	Clean sensor
No water flow	Solenoid valve defect	Replace solenoid valve
	Water supply closed	Check water supply
	Dirty sensor	Clean sensor
Equipment does not reach set temperature	Thermostat is defective	Replace thermostat
	Dirty solenoid valve	Clean solenoid valve
Equipment does not heat up	Defective circuit board	Replace circuit board
	Water level too low	Adjust water level
	Wrong temperature set	Adjust for correct temperature
	Overtemperature protection has triggered	Device from the mains supply

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