Dampfstrahler DR1

Steamer DR1

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Dear Customer

Thank you for choosing this high-quality product from BDT.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

These written instructions cannot cover all possible situations that can show up during use. For this reason, our hotline is available to answer all questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, to study the enclosed current instructions for use or to refer to our website at www.bdt-dental.de.

Manufacturer

BDT GmbH • Industriestr. 27 • 77656 Offenburg • Germany

Telephone +49 781/58 688

Fax +49 781/58 664

E-Mail info@bdt-dental.de

With reference to our quality standards

BDT ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products.

Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Instruction for use

Steamer DR-1

Item number 10013

1. Safety information

1.1 Correct usage

The steam cleaner is designed exclusively for getting rid of the wax from dentures in dental / orthodontic laboratory.

Any other use is regarded as improper use. Intended use includes the observance of all instructions in the manual.

In the final decommissioning of the BDT unit, the relevant national regulations must be observed. Questions regarding the proper disposal of the device answered BDT or dental trade.

1.2 Danger symbols and markings

In the instructions for use, the following symbols are used to indicate sources of danger to the user:

Warning This symbol indicates direct danger to the life or health of the user. Disregarding this warning may involve a serious risk to health.

Caution This symbol indicates important safety instructions. Disregarding this warning may cause damage to the unit and/or other property.

1.3 Disposal

This symbol indicates that products marked with it should not be disposed of together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through BDT or the dental trade.

2. Technical data

Mains connection	230V 50/60Hz
Rated power	2000 W
Dimensions (WxHxD)	411x500x261mm
Weight	27kg
Contents	2,75 Liter
Pressure	4,0 Bar
Protection class	I
Degree of protection	IP 33
Operable temperature range	10°C up to 40°C (50°F up to 104°F)
Use	indoor
Height above sea level	max. 1500 m ü. NN
Air humidity	10% - 80%, non-condensing

3. CE- Declaration of conformity

Through the CE marking on the rating plate we confirm compliance with the EMV Directives of the EEC 89/336/EEC, Low-voltage Directive 73/23/EEC und CE-Marking Mandatory Directive 93/68/EEC.

In the event of any modification of the product without our approval this declaration will become invalid and the warranty will extinct.

4. Commissioning of the steamer

If you should find, contrary to expectations, a transport damage to your device, please inform the company BDT or the dealer from whom you purchased it from..

Make sure before you start using the unit, that the main voltage on the rating plate (rear panel) meets your local specified operating voltage. The small steam boiler is manufactured in accordance with the technical requirements of the boiler regulation and therefore subject to permitting and monitoring free

The steam cleaner DR1 must be filled directly from the tap.

For this purpose is the 3/4 "inlet pipe on the front side of the steam cleaner should be used. Install the water inlet hose connection which has a supplied Aqua Stop included..

After connecting the inlet pipe, the sewer line needs to be connected to a sewer pipe.

Now, open the water supply (angle valve / faucet) for the steam cleaner and connect the electrical plug to the grounded outlet.

Turn on the steam cleaner at the main switch.


The water tank is filled automatically and heats up. The steam cleaner is ready in about 10 minutes.

To adjust the amount of steam the dial is provided on the solenoid valve.

Pressing the footswitch as you steam flows from the hand tool.

On the first use, air or water can escape out of the hand tool.

The water level is regulated fully automatically by means of a float switch. If the water supply fails, the power supply of the heating units turns off automatically when the water level is low.

5. Care and maintenance

The steam cleaner must be clarified every evening and it needs to be done in the following order after completion of work:

1. Switch off the electrical boiler
2. Open the on the front located ball valve and empty the entire water content.
3. Flush the steam cleaner for a half minute up to a minute with the ball valve on the main switch open.
4. Close ball valve, replenish the water until the water level is reached (red warning lamp goes off).
5. Switch off the main switch.