POLYMERISATIONSWANNE

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POLYMERIZATION BATHS

for 4 cuvettes | manually

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Instruction for use

Prüfbericht / test document

DIN EN 61010

Polymerisationswanne / Polymerization bath IVO - 4 manuell

Seriennummer / Serial number..... _____

Spannung / Voltage setting..... 230 Volt / Voltage

Verbrauchstrom / Consumption power _____ Ampere / power

Schutzleiterwiderstand 25 A Prüfstrom /
Protective conductor resistance 25A testing current > _____ Ohm

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Dear Customer

Thank you for choosing this high-quality product from BDT.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

These written instructions cannot cover all possible situations that can show up during use. For this reason, our hotline is available to answer all questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, to study the enclosed current instructions for use or to refer to our website at www.bdt-dental.de.

Manufacturer

BDT GmbH • Sonnhalde 1 • 77749 Hohberg • Germany

Telephone +49 7808 914996-0

Fax +49 7808 914996-99

E-Mail info@bdt-dental.de

With reference to our quality standards

BDT ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products.

Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Instruction for use

Polymerization baths

Artikelnr. 10006

1. Safety information

1. Correct usage

The polymerization bath is intended only for polymerizing plastic prostheses in the dental laboratory.

Any other use shall be regarded as not intended. Intended use shall also include compliance with any and all information contained in the instructions for use.

The corresponding regulations of the respective country must be complied with once the BDT equipment is finally taken out of service. BDT or the specialized dental retailer can answer any questions with regard to proper disposal of the equipment.

2. Danger symbols and markings

In the instructions for use, the following symbols are used to indicate sources of danger to the user:

Warning This symbol indicates direct danger to the life or health of the user. Disregarding this warning may involve a serious risk to health.

Caution This symbol indicates important safety instructions. Disregarding this warning may cause damage to the unit and/or other property.

3. Disposal

This symbol indicates that products marked with it should not be disposed of together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through BDT or the dental trade.

2. Description and usage

The polymerization bath is used for polymerizing plastic prostheses in the dental laboratory.

The polymerization bath consists of a stainless steel bath with electronic water level control. The temperature can be continuously controlled. As a result of indirect heating there is no furring-up of heating elements.

NOTE: Please note that optimal polymerization can only be ensured with regularly serviced equipment. In this case national regulations must also be taken into consideration.

3. Technical data

Mains connection	230V	50/60Hz
Rated power	2500W	
Degree of protection	IP 33	
Discharge hose	Ø22mm	
Dimensions (WxHxD)	530x300x470mm	
Weight	14kg	
Operable temperature range	10°C up to 40°C (50°F up to 104°F)	
Use	indoor	
Height above sea level	max. 1500 m ü. NN	
Air humidity	10% - 80%, non-condensing	

4. Commissioning of the suction unit

Immediately notify BDT or the specialized dealer where you purchased the equipment should you, contrary to expectations, determine any transport damage to your equipment.

Prior to commissioning make certain that the mains voltage agrees with the operating voltage indicated on the model identification plate (back of the equipment).

After setting up the equipment the waste water hose must be attached to the waste system). After the equipment is connected to the electrical mains (230 V) the standpipe is inserted into the overflow valve. The equipment is switched on by pressing the main switch. Now fill in water in the pool up to the end of the standpipe. Heating is possible only when the prescribed water level has been reached. The desired temperature – which must be **below the boiling point** – is now adjusted on the thermostat. To avoid the development of steam, the Floating ball can be used.

(cf. Fig. 1)

For safety reasons the heating switches off automatically if the water level is too low.


The water intake must be turned off each evening and the equipment disconnected from the electrical mains.

Do not fill in combustible liquids (fire hazard).

picture 2

5. Care and maintenance

If necessary wipe off the housing with a dry or dampened cloth. Do not use any strong solvents.

The stainless steel bath should preferably be cleaned with a basic stainless steel cleaner and Breglit stainless steel care spray.

The water level sensor must be cleaned once per month.

Repairs should only be carried out by authorized technical personnel.

Only open the housing if the equipment has been disconnected from the electricity mains.

6. CE- Declaration of conformity


Through the CE marking on the rating plate we confirm compliance with the EMV Directives of the EEC 89/336/EEC, Low-voltage Directive 73/23/EEC und CE-Marking Mandatory Directive 93/68/EEC.

In the event of any modification of the product without our approval this declaration will become invalid and the warranty will extinct.

7. Error search

Problem	Possible Cause	Remedy
Green control lamp not lit	Mains cable unplugged	Plug in mains cable
	Fuse defect	Replace fuse
	No mains voltage	Check mains voltage
Water always flows in	Solenoid valve does not close	Clean solenoid valve
	Standpipe is improperly or not installed	Check proper fit of standpipe
	Sensor set at wrong height	Set correct sensor height
	Dirty sensor	Clean sensor
No water flow	Solenoid valve defect	Replace solenoid valve
	Water supply closed	Check water supply
	Dirty sensor	Clean sensor
Equipment does not reach set temperature	Thermostat is defective	Replace thermostat
	Dirty solenoid valve	Clean solenoid valve
Equipment does not heat up	Defective circuit board	Replace circuit board
	Water level too low	Adjust water level
	Wrong temperature set	Adjust for correct temperature

IVO-4 manuell 230 V

- 20018 Standrohr 28x115mm
pipe 28x115mm
- 20019 Standrohrventil
drainage
- 50031 Dichtung 33x43x3mm
sealing 33x43x3
- 70013 Mutter R3/8
nut R3/8
- 50051 Verschlusschraube
stopper R3/8
- 20180 Einlegeblech 476x275mm
Inlay sheet 476x275mm
- 20051 Becken mit Heizung 230 V /2500W
basin with basin heating 230V 2500W
- 20190 Isolierplatte 505x295mm
isolation 505x295mm
- 60081 Dichtung 16x22x2 mm
sealing 16x22x2mm
- 60052 Einlaufwinkel R3/8xR3/8
pipe 90° R3/8xR3/8
- 70015 Zahnscheibe
star washer
- 30024 Erdungsschild
ground shield
- 70057 Schraube M4x8mm
screw M4x8mm
- 60017 Mutter R1
nut R1
- 60010 Schlauchklemme
hose bandage
- 20081 Fühlerwasserstand
water level feeler
- 70094 Schraube M3x8mm
screw M3x8mm
- 20010 Ablaufschlauch
waste water hose
- 70023 Schraube 2x10mm
screw 2x10mm
- 70002 Lüftungsgitter
ventilation screen
- 30016 Verschraubung-Zuleitung
no pull fastner
- 30077 Verschlussfülle
stopper Z1
- 30011 Schelle
no pull fastener for hoses
- 30012 Netzanschlußstecker mit Kabel
plug and fead
- 20169 Boden
cover
- 20022 Gerätefüße
housing supports
- 70083 Schrauben M4x12mm
screw M4 x12mm
- 20030 Gehäuse
housing
- 40042 Platine
bord
- 30010 Hauptschalter
main power switch
- 30017 Sicherungshalter
fuse holder
- 30002 Sicherung 5x20x12,5A
fuse 5 x20 x12,5A
- 30007 Thermostat
thermostate with temperature sensor